






























# restor3d

## Symbols Glossary

Symbol	Definition	Reference / Standard
	Batch Code	ISO 15223-1:2021 5.1.5
	Catalog number	ISO 15223-1:2021 5.1.6
	Serial Number	ISO 15223-1:2021 5.1.7
	Patient Identification	ISO 15223-1:2021 5.7.3
	Health Care Center (name and address) or Doctor	ISO 15223-1:2021 5.7.5
	Quantity	N/A
	Do not re-use	ISO 15223-1:2021 5.4.2
	Consult Instructions For Use.	ISO 15223-1:2021 5.4.3
	Consult Instructions For Use. Available on the restor3d website at <a href="http://www.restor3d.com">www.restor3d.com</a>	ISO 15223-1:2021 5.4.3
	Temperature limitation	ISO 15223-1:2021 5.3.7
	Keep dry	ISO 15223-1:2021 5.3.4
	Keep away from sunlight	ISO 15223-1:2021 5.3.2

	Date of manufacture	ISO 15223-1:2021 5.1.3
	Manufacturer	ISO 15223-1:2021 5.1.1
	Use by date, use by, expiry date, expiration date	ISO 15223-1:2021 5.1.4
	Sterilized using ethylene oxide	ISO 15223-1:2021 5.2.3
	Sterilized using irradiation	ISO 15223-1:2021 5.2.4
	Sterilized using vaporized hydrogen peroxide	ISO 15223-1:2021 5.10
	Sterilized using vaporized hydrogen peroxide	N/A – Legacy VHP symbol
	Non-Sterile	ISO 15223-1:2021 5.2.7
	Do not re-sterilize	ISO 15223-1:2021 5.2.6
	MR Unsafe	ASTM F2503 – 23
	MR Conditional	ASTM F2503 – 23
<b>Rx Only</b>	For prescription use only	21 CFR 801.15(c)(1)(i)F 21 CFR 801.109
	Do not use if the package has been damaged or opened and consult instructions for use	ISO 15223-1:2021 5.2.8
	Single sterile barrier System	ISO 15223-1:2021 5.2.11

	Double sterile barrier System	ISO 15223-1:2021 5.2.12
	Single sterile barrier system with protective packaging inside	ISO 15223-1:2021 5.2.13
	Single sterile barrier system inside protective packaging	ISO 15223-1:2021 5.2.14
	Caution (when operating the device)	ISO 15223-1:2021 5.4.4